Creatinine is the product of phosphocreatine decomposition. It is produced endogenously and released into body fluids at a constant rate. Its plasma concentration is maintained within narrow limits, predominantly by glomerular filtration. Consequently, both the plasma creatinine concentration and its renal clearance (creatinine clearance) have been used as markers of glomerular filtration rate.

PERFORMANCE CHARACTERISTICS

	Operating range	
	Serum	0.28 a 12.00 mg/dL
	Urine	3.15 a 321.35 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

	Sensitivity	
	Detection Limit	Quantification Limit
	Serum	0.11 mg/dL 0.28 mg/dL
Urine	0.018 mg/dL	3.15 mg/dL

	Analytical Specificity		
	Hemoglobin	Bilirubin	Triglycerides
Serum	500 mg/dL	10 mg/dL	500 mg/dL
Urine	180 mg/dL	6 mg/dL	-

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum	Urine
Number of Samples	40 in duplicate	20 in duplicate
Regression Equation	y = 0.998x + 0.040	y = 1.002x - 2.343
Correlation Coefficient (R)	0.9958	0.9995

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 1.0 mg/dL and 4.0 mg/dL was 3.8% and 0.8%, respectively. Urine: by applying the regression equation, the total systematic error estimated for decision levels of 1500 mg/24h and 2000 mg/24h was 0.04% e 0,08%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Serum Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
0.96	80	0.004	0.40	0.005	0.50
4.07	80	0.012	0.30	0.017	0.40
10.41	80	0.125	1.20	0.132	1.30

Urine Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
93.69	80	1.699	1.80	2.202	2.40
182.10	80	3.921	2.20	4.176	2.30

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

	Men		Women	
	mg/dL	µmol/L	mg/dL	µmol/L
	0,9 - 1,3	80 - 115	0,6 - 1,1	53 - 97
Urine	mg/kg/day	µmol/kg/day	mg/kg/day	µmol/kg/day
	14 - 26	124 - 230	11 - 20	97 - 177

Clearance (mL/min) - Adults	52,1 – 110
-----------------------------	------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µmol/L):

Creatinine (mg/dL) x 88,4 = Creatinine (µmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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